

Product Information: Similac® Nipples

For more information, contact your Abbott Nutrition Representative or visit
www.abbottnutrition.com

- Provides the flexibility to choose a ready-for-use nipple and ring that match each infant's sucking strength or need.
- **Infant** (latex-free)—for infants with developed sucking strength.
- **Premature** (latex-free)—for preterm infants with less developed sucking strength.
- **Orthodontic** (contains latex)—for newborn infants, to stimulate oral exercise while feeding.
- **Slow Flow** (latex-free)—single hole for a more controlled flow rate during feeding.
- **Cleft Palate Assembly** (latex-free)—uniquely designed to meet the feeding needs of the infant with cleft palate. Requires minimal sucking. Formula flow can be regulated by gently squeezing the side walls of the attached bottle.
- Nipples and rings individually prepackaged and ready for use.
- Nipples and rings form a seal to prevent leakage during feeding.
- Disposable—not for reuse.
- Easy to use.

Safety Precautions

- DO NOT use as a pacifier, as it presents a choking hazard.
- DO NOT reuse, as reuse can cause deterioration or separation of the rubber, presenting a choking hazard.
- Keep out of reach of children.

Availability

List Number	Item
00079	Similac Infant Nipple & Ring / 250 ct
00094	Similac Premature Nipple & Ring / 250 ct
50512	Similac Orthodontic Nipple & Ring / 250 ct
53894	Similac Slow Flow Nipple and Ring / 250 ct
56054	Cleft Palate Assembly / 3 ct
56366	Similac Slow Flow Nipple & Ring Kit / 50 ct

For more information, contact your Abbott Nutrition Representative or visit
www.abbottnutrition.com