# RECOMMENDED STORAGE GUIDELINES FOR ABBOTT NUTRITION PRODUCTS

# RECOMMENDED STORAGE TEMPERATURES FOR SEALED CONTAINERS†

### Ready-to-use or concentrated liquid

- Most desirable: between 13°C and 24°C.
- Acceptable: between 0°C and 35°C.
- Do not freeze.

Prolonged exposure to temperatures below 0°C or to direct heat above 35°C could affect the physical consistency of the product and consequently could also affect its appearance, flavour, and other sensory attributes. Product is sterile and should remain of sound nutritional quality.



#### Powder

- Store in dry, cool area.
- Do not refrigerate.

# GUIDELINES FOR PREPARING AND STORING FORMULA

#### **Human Milk Fortifiers**

See individual product pages for storage information on Similac® Human Milk Fortifier Special Care® 30 (with iron), Similac® Human Milk Fortifier, and LiquiProtein™.

# Ready to use

<u>Cans and bottles</u>: after opening, cover, refrigerate, and use within 48 hours.

945-mL reclosable plastic bottle: re-cap, refrigerate, and use within 72 hours.

## Concentrated Liquid

- Open can (undiluted): cover, refrigerate, and use within 48 hours.
- Prepare as per label instructions.
- Pour into individual feeding bottles or cups, cover, and store in refrigerator.
- Use prepared formula within 48 hours.

#### Powder

- Prepare as per label instructions.
- Pour into individual feeding bottles or cups, cover, and store in refrigerator.
- · Use within 24 hours.
- · Cover opened container after each use.
- Use product within 1 month of opening container.

If ready-to-use, undiluted or prepared formula is left at room temperature, it should be used within 2 hours.<sup>1</sup>

# RECOMMENDED STORAGE GUIDELINES FOR ABBOTT **NUTRITION PRODUCTS**

# GUIDELINES FOR PREPARING AND STORING FORMULA (continued)

If product is consumed directly or by straw from can or bottle:‡

### Infant formula

- · Dispose of any product remaining in the bottle within 1 hour from beginning of feed (first time baby's mouth comes into contact with nipple).
- Do not reuse 59-mL (2-oz) plastic bottles.

# PediaSure® products and Similac® Mom (adult nutritional product)

- · Use or refrigerate within 1 hour.
- If refrigerated, use product within 24 hours.

# HANG TIME FOR ABBOTT NUTRITION INFANT AND PEDIATRIC PRODUCTS

#### Ready to use (commercially sterile)

• 4-8 hours when carefully poured from bottle/can into an enteral-feeding setup.2‡

## Reconstituted products

Includes: Concentrated liquid, powder, any formula to which a modular component has been added, or human milk that has been fortified.

- · Prepare using aseptic technique.
- Hang for no more than 4 hours.<sup>2</sup>

Each health care facility is responsible for developing/ following procedures for using safe techniques in the preparation and use of tube feedings and in setting quality assurance protocols for monitoring those techniques.

# EXPIRATION DATES

- To assure the highest quality, use the formula by the date shown on the product.
- The vitamin content and the physical stability of the product cannot be guaranteed beyond the expiry date because both may degrade with time.

# WARNINGS

- Never use a microwave oven to warm formula: serious burns can result.
- · Abbott Nutrition products are not intended for parenteral use.

Tints produce exposure or programment of the progr

<sup>†</sup> For highest quality product.

<sup>‡</sup> This practice exposes the product to significant amounts of bacteria